




## Imaging Chamber

## Imaging Chamber CG

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
## Application technologies

Laser Scanning Confocal Microscopy (LSM), Total Internal Reflection Fluorescence (TIRF), Differential Interference Contrast (DIC), Fluorescence Correlation Spectroscopy (FCS), Fluorescence Resonance Energy Transfer (FRET), Fluorescence Recovery After Photobleaching (FRAP), Low intensity fluorescence, Fluorescence In-Situ Hybridization (FISH), Immunohistology

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Konfokale Laser Scanning Mikroskopie (LSM), Totale Internal Reflection Fluorescence (TIRF), Differentieller Interferenz Kontrast (DIC/Nomarski), Fluoreszenz Korrelationsspektroskopie (FCS), Fluorescence Resonance Energy Transfer (FRET), Fluorescence Recovery After Photobleaching (FRAP), Schwach fluoreszierende Farbstoffe/Proteine, Fluoreszenz In-Situ Hybridisierung (FISH), Immunhistologie

**Imaging Chamber**  
**Imaging Chamber CG**

## Features & benefits

Imaging Chambers are available with microscope slides or cover glass slides (Imaging Chamber CG) as bottom structure. Chambers and lids are made from high performance polymers and are sealed against the glass bottoms by biological inert silicone adhesive. At the end of the experiments the glasses can be detached from the polymer chambers with minimal glue residue.

The high performance polymer of chamber bodies and lids withstands temperatures up to 90°C and is resistant against acetone. Therefore the full products are suitable for in-situ hybridization technologies as well as for a wide range of fixation and staining protocols.

Imaging Chambers CG are designed for high resolution live cell microscopy techniques. Imaging Chambers with microscope slides can be used with long distance objectives. After disassembly the slides can be covered with cover glasses and are then also suitable for high resolution immersion microscopy.

Cell adhesion, spreading and distribution on top of the glass surface is improved due to our proprietary plasma surface modifications.

Imaging Chambers are available in 1 well, 2 well, 4 well and 8 well format.

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

## Product highlights

- **extraordinary planicity** of the imaging plane, even during temperature shifts
- high temperature stability, **suitable for FISH**
- high chemical resistance, **fixation with acetone** enabled
- **bottom glasses** can be obtained **with minimal glue residues**

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
Thickness microscope slides	1 mm
Dimensions microscope slides	26 x 76 mm
Microscope slides: soda lime glass,	hydrolytic class 3
Thickness cover glass	170 µm
Dimensions cover glass	26 x 58 mm
Cover glasses: borosilicate glass, hydrolytic class 1	
Surface treatment	Tissue Culture (TC)
Refractive index	1.52
Abbe's number	55
Planicity/Flatness	≤ 5µm


## Imaging Chamber Imaging Chamber CG

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### General information


All Imaging Chambers are delivered with lid and are placed on a Slide Tray Plate.

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### Imaging Chamber CG adapter

Re-usable metal clip to extend Imaging Chamber CG to full length of microscope slides (76 mm). Usefull if your microscope stage insert is prepared to carry standard slides (26 mm x 76 mm).

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**Imaging Chamber**  
**Imaging Chamber CG**

### Order information / Bestellinformationen

Art.No.:	Imaging Chamber	Qty/box	Art.No.:	Imaging Chamber CG	Qty/box
8011-4	1 well - 4/bag	4	8001-4	1 well - 4/bag	4
8011-16	1 well - 4/bag	16	8001-16	1 well - 4/bag	16
8011-80	1 well - 4/bag	80	8001-80	1 well - 4/bag	80
8012-4	2 well - 4/bag	4	8002-4	2 well - 4/bag	4
8012-16	2 well - 4/bag	16	8002-16	2 well - 4/bag	16
8012-80	2 well - 4/bag	80	8002-80	2 well - 4/bag	80
8014-4	4 well - 4/bag	4	8004-4	4 well - 4/bag	4
8014-16	4 well - 4/bag	16	8004-16	4 well - 4/bag	16
8014-80	4 well - 4/bag	80	8004-80	4 well - 4/bag	80
8018-4	8 well - 4/bag	4	8008-4	8 well - 4/bag	4
8018-16	8 well - 4/bag	16	8008-16	8 well - 4/bag	16
8018-80	8 well - 4/bag	80	8008-80	8 well - 4/bag	80
8100-4	DIC Lid	4			
8200-1	Imaging Chamber CG adapter	1			

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

# Imaging Chamber / Imaging Chamber CG

Instructions for handling and chamber  
disassembly

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

All Imaging Chamber and Imaging Chamber CG products are designed to allow disassembly of the chamber body and the glass bottom at the end of experiments / investigations. Due to their physical characteristics different strategies for microscope slide bottoms and the cover glass bottoms have to be considered. Also the chamber design (one well versus multiwell) determines the optimal technique for disassembly.

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
## Handling / Handhabung

During transport we recommend to hold the chambers where they are supported by a wall or to grip the microscopic slide (Imaging Chamber products).

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## Disassembly / Demontage Imaging Chamber 1 und Imaging Chamber 1 CG


To disassemble the glass bottom from Imaging Chamber 1 or Imaging Chamber 1 CG you only need to compress the two sidewalls of the chamber. The shear force detaches the glass from the chamber.

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## Imaging Chamber / Imaging Chamber CG

Instructions for handling and chamber  
disassembly

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



To disassemble the microscope slide from Imaging Chamber 2, 4 or 8 you have to hold the slide at the outer frosted end close to the chamber while lifting the chamber beginning from this end with the other hand.

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## Disassembly / Demontage Imaging Chamber 2, 4 und 8 CG Cover Glass / Deckglas

To disassemble the cover glass from Imaging Chamber 2, 4 or 8 CG you will need a tool like a spatula. Place the imaging Chamber on a clean and flat ground. Lift the Chamber approximately 1 mm at one end of the Chamber from the ground. Gently press the cover glass with the spatula away from the chamber body close to the inner Chamber walls. Proceed well by well until the glass comes off.

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

## Assembly

(Patent Pending 10-2011-001328, SPL)


## 1. Eco-friendly and more convenient design

- 1) Use of no adhesives and adaption of mechanical assembly help you detach the chamber off the slide as easily ever without breaking it.
- 2) Gives you cyto-toxic freedom by non use of adhesives and you don't have to be afraid of destroying your cultured output while detaching the chamber off the slide

## 2. Excellent test results and secure sealing

- 1) Confocal Image: you can get the crystal clear image without background (tested by C university, domestic)


Phase contrast

DAPI

RITC

Merge

- 2) Strong sealing free of leakage

## 3. Specification of Product


Type	Cat. No.	Chamber type Number / Color	Material Chamber / Slide / Holder	Vol. / Well	Sterilized	Units per Sleeve / Case
	30104	4 / Clear	PS / Glass / PP	0.5 ~ 1.3 ml	+	6 / 12
	30114	4 / Black	PS / Glass / PP	0.5 ~ 1.3 ml	+	6 / 12
	30124	4 / White	PS / Glass / PP	0.5 ~ 1.3 ml	+	6 / 12
	30108	8 / Clear	PS / Glass / PP	0.2 ~ 0.6 ml	+	6 / 12
	30118	8 / Black	PS / Glass / PP	0.2 ~ 0.6 ml	+	6 / 12
	30128	8 / White	PS / Glass / PP	0.2 ~ 0.6 ml	+	6 / 12

## Application technologies

Laser Scanning Confocal Microscopy (LSM), Total Internal Reflection Fluorescence (TIRF), Differential Interference Contrast (DIC), Fluorescence Correlation Spectroscopy (FCS), Fluorescence Resonance Energy Transfer (FRET), Fluorescence Recovery After Photobleaching (FRAP), Low intensity fluorescence, Microinjection / Micromanipulation, Cytometry

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Konfokale Laser Scanning Mikroskopie (LSM), Totale Internal Reflection Fluorescence (TIRF), Differentieller Interferenz Kontrast (DIC/Nomarski), Fluoreszenz Korrelationsspektroskopie (FCS), Fluorescence Resonance Energy Transfer (FRET), Fluorescence Recovery After Photobleaching (FRAP), Schwach fluoreszierende Farbstoffe/Proteine, Mikroinjektion / Mikromanipulation, Zytometrie

**Imaging Dish CG**  
**Imaging Dish Stage Frame**

## Features & benefits

### 35 mm petri dish with cover glass bottom

Imaging Dish CG are designed for high resolution microscopy techniques. The 35 mm petri dishes provide a central cover glass bottom with a diameter of 18 mm which supports direct inverse microscopy of living and fixed cells. The robust injection moulded dish body allows an extraordinary planicity over the entire imaging area.

Cell adhesion, spreading and distribution on top of the glass surface is improved due to our proprietary plasma surface modifications.

The cover glass area and cell cultivation zone is kept 2 mm below the basement level of the petri dish body. Therefore a small volume reservoir is created which allows the concentration of the cells on top of the glass area during seeding and reduction of antibodies or dyes during staining applications.

Special attention was given to the polygonal grip zone of the dishes for improved handling and the prominent macroscopic orientation marks (N-O-S-W) on the dish bottom. The combination of the polygonal grip zone of the dishes with our Imaging Dish Stage Frame allows exact repositioning of the dishes for repeated imaging sessions on the microscope.

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
## Product highlights


- **extraordinary planicity** of the imaging plane, even during temperature shifts
- easy handling due to **ergonomic grip zone and lid design**
- Imaging Dish  $\mu$ Grid with **microscopic grid for cytometry** and recovery of ROI (Region Of Interest)
- **Imaging Dish Stage Frame** for exact and easy positioning and orientation on the microscope stage

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

Dish diameter	35mm
Imaging area / diameter	18mm
Total volume	7ml
Suggested cell seeding volume	750µl
Suggested working volume	2ml
Distance dish bottom / focus plane	400µm
Surface treatment	Tissue Culture (TC)
Lid	yes
Sterile	yes
Planicity/Flatness	≤ 5µm


## Imaging Dish CG Imaging Dish Stage Frame


### Imaging Dish µGrid

Central high quality micro grid, grid lines 9 mm long, 3 µm wide and 0.3 µm deep, grid aligned to macroscopic orientation marks  
 Inner grid line distance: 300 µm  
 Outer grid line distance: 600 µm  
 Easy detection of grid in brightfield and phase contrast, no impact on fluorescence image acquisition

### Imaging Dish Stage Frame

ID Stage Frame single, 52 x 76 mm, usable on stages with microscope slide holders. ID Stage frame double, 124,76 x 85,48 mm, usable on stages with multiwell inserts. ID stage frames are made from aluminium. Please clean with soft tissues if necessary to avoid scratches.

### General information

Cover glasses made from borosilicate glass, hydrolytic class I with low intrinsic fluorescence. Glasses with 145 µm thickness (1.0) and 170 µm (1.5) available. Dish body and lid are made from polystyrene (PS). Do not use chemicals not tolerated by PS like acetone. Imaging dishes can be used between -20° and 50°C.

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### Tip for DIC/Nomarski

Superior Results for DIC Imaging can be achieved when re-usable glass cover lids from Bioswisstec are used in combination with Imaging Dishes.

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## Application technologies

Laser Scanning Confocal Microscopy (LSM), Total Internal Reflection Fluorescence (TIRF), Differential Interference Contrast (DIC), Fluorescence Correlation Spectroscopy (FCS), Fluorescence Resonance Energy Transfer (FRET), Fluorescence Recovery After Photobleaching (FRAP), Low intensity fluorescence, Fluorescence In-Situ Hybridization (FISH), Immunohistology, High Content Analysis/Screening (HCA/HCS)

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Konfokale Laser Scanning Mikroskopie (LSM), Total Internal Reflection Fluorescence (TIRF), Differentieller Interferenz Kontrast (DIC/Nomarski), Fluoreszenz Korrelationsspektroskopie (FCS), Fluoreszenz Resonance Energy Transfer (FRET), Fluoreszenz Recovery After Photobleaching (FRAP), Schwach fluoreszierende Farbstoffe/Proteine, Fluoreszenz In-Situ Hybridisierung (FISH), Immunhistologie, High Content Analysis/Screening (HCA/HCS)

**Imaging Plate CG**  
**24 well / 96 well**

## Features & benefits

Imaging Plates CG are black multiwell plates compliant to the SBS (Society for Biomolecular Screening) standard with a bottom made from borosilicate cover glass. The plate bodies are made from polystyrene and are tightly sealed against the glass bottoms by adhesive.

Imaging Plates CG can be used for high resolution live cell microscopy. The plates are also suitable for a wide range of fixation and staining protocols (formaldehyde, glutaraldehyde, alcohols, acetone). Acetone should not be applied for periods longer than 10 minutes since the plate body is made from polystyrene and the plate body can detach from the glass bottom.

Imaging Plates CG can be used at temperatures between -20°C and 50°C.

Cell adhesion, spreading and distribution are improved due to our well established plasma surface modification technologies.

The cover glass plate bottoms make the plates suitable for the application of high resolution immersion objectives. Water, glycerine and oil can be used as immersion medium.

Imaging Plates CG are available in 24 well and 96 well format.

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
## Product highlights


- **extraordinary planicity** of the imaging plane
- **high resolution** imaging for **throughput applications**
- 96 well plate with "**low-skirt design**", full bottom area can be used with immersion objectives
- **Advanced Tissue Culture** surface quality for cell adhesion and spreading

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
### Imaging Plate CG 24 well / 96 well

Article number / Artikelnummer	5231	5241
Format	24 well	96 well
Plate length	127,76 mm	127,76 mm
Plate width	85,48 mm	85,48 mm
Inner well diameter (imaging plane)	13.2 mm	6 mm
Distance well to well center	18 mm	9 mm
Distance center well A1 plate border (long side)	15.74 mm	11.24 mm
Distance center well A1 plate border (small side)	18.88 mm	14.38 mm
Distance plate bottom to image focus plane	2.5 mm	400 µm
Inner well flatness	≤ 10 µm	≤ 10 µm
Well to well flatness	≤ 50 µm	≤ 50 µm
Total well volume	1880 µl	428 µl
Suggested working volume	500-1000 µl	100-200 µl

The polymer bodies of all labware react to temperature shifts with changes of their dimension. Please ensure that you can achieve equilibrated environmental temperature conditions during microscopy to avoid drift of the focal plane.

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#### General information

All Imaging Plates are shipped single packed with lid.  
All lids provide condensation rings and spacers for optimal ventilation.

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
## Application technologies

- Horizontal handling of microscope slides (e.g. for Comet-Assays)
- Horizontal handling of microarray slides
- Staining / storage / incubation of slides
- Transport and handling of chambered slides and slide flasks

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**Slide Tray Plate PS / PS TC**  
**Slide Tray Plate PP**

## Features & benefits

### Four well plate in SBS footprint format

Up to four slides can be handled in parallel. Each well is marked (A-D) and contains an access zone for easy slide removal.

Lids and plates are optimized for stacking and a grip zone makes the handling ergonomic.

Two material variants are offered:  
Transparent PS (polystyrene) for cell culture incubation or Chemical resistant PP (polypropylene).

Slide Tray Plate PS TC have a surface treatment (amine groups) for optimal cell adhesion and distribution and are suited for direct cell cultivation of adherend cells.

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## Order information / Bestellinformationen

Art. No.	Product	Surface	Qty/Bag	Qty/Box
1145-40	Slide Tray Plate PS, sterile	hydrophobic	10	40
1245-40	Slide Tray Plate PS TC, sterile	Tissue Culture	10	40
1345-40	Slide Tray Plate PP	hydrophobic	10	40

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## Application technologies

UV Spectroscopy, Transmission and detection of Optical Density (OD) between 220 nm and 280nm, Evaluation of DNA and RNA purity, label free quantitation of DNA, RNA and small molecules.

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UV Plate

## Features & benefits

### UV spectroscopy in 96 well format and throughput

Due to its ultrathin film bottom (25 µm) made from a high performance polymer this microwell plate shows excellent light transmission light in the UV A and UV B range.

In combination with an UV plate reader the purity of DNA and RNA can be detected. Furthermore the concentration of DNA and RNA can be detected without addition of any stains and markers.

The plate is designed in accordance with SBS guidelines and is suitable for automated throughput screenings.

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
## Technical data / Technische Daten

Format	96 well
Inner well diameter	6mm
Total volume per well	428 µl
Suggested working volume	50 - 350 µl
Surface treatment	hydrophilic
Lid	no
Sterile	no

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prices in CHF excl. VAT 8% ex works Schaffhausen, supply with the most economic transport system

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Bioswisstec AG